# Non-invasive Blood Pressure Accessories



## Non-invasive Blood Pressure Tubing

- Not deform even if for long-term use
  Reasonable hardness avoids severe bending, ensures the safety of deflation
  Reasonable hardness ensures well transmission of blood pressure signal
  TPU ensures good air tightness and long life
  Latex free, PVC free
  Good blocompatibility, free from biological hazard to skin

For BeneVision, BeneView, iPM, iMEC, uMEC series monitors, VS-600/900, BeneHeart defibrillator				
Picture	Part No.	Description	Purchasing Unit	
	6200-30-09688	NIBP Tubing, Adu/Ped/Inf, with connectors, 3 m	Each	
	115-012522-00	NIBP Tubing, Adu/Ped/Inf, with connectors, 3 m (For iMEC, uMEC series patient monitor)	Each	



6200-30-11560 NIBP Tubing, Neo, with connectors, 3 m Each

#### For VS-800, MEC, PM series monitors



Part No.	Description	Purchasing Unit
509B-30-06259	NIBP tubing, Adu/Ped/Inf, 3 m	Each



509B-30-06260 NIBP tubing, Neo, 3 m Each

## Adapted with the tubing (6200-30-09688 and 509B-30-06259) and the neonatal cuffs



Picture	Part No.	Description	Purchasing Unit
)28	040-000688-00	NIBP Cuff Tubing Adapter (Adult tubing to Neonate cuff)	Each

# Non-invasive Blood Pressure Cuff (Reusable)

#### CM1200 Series

- Soft and comfortable. Low hazard to skin even if a long-term use
  Easy to clean. The cuff wrap can not be damped or stained by liquid if duly cleaned
  Pilling-proof. Not deform even if for long-term use
  TPU bladder ensures good air tightness and long life
  Latex free, PVC free
  Good biocompatibility, free from biological hazard to skin



## Connected with the tubing 6200-30-09688, 509B-30-06259 and 115-012522-00

.ormected w	ntil the tubing 0200-5	0=09000, 3090=30=00239 and 113=012322=00	
Model	Part No.	Description	Purchasing Unit
CM1200	115-002480-00	Reusable cuff, Small Inf, 7-13 cm Connected with the tubing 5098-30-06260 and 6200-30-11560	Each
CM1201	0010-30-12157	Reusable cuff, Inf, 10-19 cm, with connector	Each
CM1202	0010-30-12158	Reusable cuff, Child, 18-26 cm, with connector	Each
CM1203	0010-30-12159	Reusable cuff, Adu, 25-35 cm, with connector	Each
CM1204	0010-30-12160	Reusable cuff, Large Adu, 33~47 cm, with connector	Each
CM1205	0010-30-12161	Reusable cuff, Thigh, 46-66 cm, with connector	Each

## CM1300 Series

- Easy to clean. No bladder, can be cleaned and disinfected directly
- Color-coded by size for different patient
- Avoid inaccurate measurement. No-bladder design avoids inaccurate measurements due to bladder twisting

- messurements due to alaboret (wissing Pillling-proof. Not deform even if for long-term use Round angle design improves safety Easy to use. The index line for proper size and placement Latex free, PVC free
- Good biocompatibility



### Connected with the tubing 6200-30-09688, 509B-30-06259 and 115-012522-00 $\,$

Model	Part No.	Description	Purchasing Unit
CM1300	040-000968-00	Reusable bladderless cuff, Small Inf, 7-13 cm Connected with the tubing 5098-30-06260 and 6200-30-11560	Each
CM1301	040-000973-00	Reusable bladderless cuff, Inf, 10-19 cm, with connector	Each
CM1302	040-000978-00	Reusable bladderless cuff, Child, 18-26 cm, with connector	Each
CM1303	040-000983-00	Reusable bladderless cuff, Adu, 25-35 cm, with connector	Each
CM1304	040-000988-00	Reusable bladderless cuff, Large Adu, 33–47 cm, with connector	Each
CM1305	040-000993-00	Reusable bladderless cuff, Thigh, 46–66 cm, with connector	Each
CM1306	115-015930-00	Reusable bladderless cuff, Adu, long wrap, 25-35 cm, with connector	Each
CM1307	115-015931-00	Reusable bladderless cuff, Large Adu, long wrap, 33-47 cm, with connector	Each

# Non-invasive Blood Pressure Cuff (Single Use) CM1500 Series (Adult and Child)

- Soft and comfortable. Low hazard to skin even if a long-term use
   Avoid inaccurate measurement. No bladder design avoids inaccurate measurements due to bladder twisting

- For single patient only, avoiding cross infection
   Round angle design improves safety
   Easy to use. The index line for proper size and placement
- Latex free, PVC free
- Good biocompatibility, free from biological hazard to skin



### Connected with the tubing 6200-30-09688, 509B-30-06259 and 115-012522-00

Connected with the tabing 0200-30-05000, 305b-30-00235 and 113-012322-00				
Model	Part No.	Description	Purchasing Unit	
CM1501	001B-30-70697	Single-patient use cuff, Inf, 10-19 cm	10 pcs/box	
CM1502	001B-30-70698	Single-patient use cuff, Child, 18-26 cm	10 pcs/box	
CM1503	001B-30-70699	Single-patient use cuff, Adu, 25-35 cm	10 pcs/box	
CM1504	001B-30-70700	Single-patient use cuff, Large Adu, 33-47 cm	10 pcs/box	
CM1505	001B-30-70701	Single-patient use cuff, Thigh, 46-66 cm	10 pcs/box	
CM1506	115-016969-00	Single-patient use cuff, Adu, long wrap, 25-35 cm	10 pcs/box	
CM1507	115-016970-00	Single-patient use cuff, Large Adu, long wrap, 33-47 cm	10 pcs/box	

## CM1500 Series (Neonatal)

- Soft and comfortable. Low hazard to skin even if a long-term use
  Easy to clean. The stain can be cleaned by wet cloth
  Avoid inaccurate measurement
  No-bladder design avoids inaccurate measurements due to bladder twisting
- For single patient only, avoiding cross infection

  Transparent design for observing the patient skin condition conveniently

  Easy to use. The index line for proper size and placement

- Latex free, PVC free
   Good biocompatibility, free from biological hazard to skin



## Connected with the tubing 6200-30-11560 and 509B-30-06260

Model	Part No.	Description	Purchasing Unit
CM1500A	0018-30-70692	Single-patient use cuff, Neo, 3.1-5.7 cm	20 pcs/box
CM1500B	001B-30-70693	Single-patient use cuff, Neo, 4.3-8.0 cm	20 pcs/box
CM1500C	001B-30-70694	Single-patient use cuff, Neo, 5.8-10.9 cm	20 pcs/box
CM1500D	001B-30-70695	Single-patient use cuff, Neo, 7.1-13.1 cm	20 pcs/box
CM1500E	001B-30-70696	Single-patient use cuff, Neo, 8-15 cm	20 pcs/box

Picture	Part No.	Description	Purchasing Unit
	040-002379-00	NIBP tubing	Each
Juna	115-020374-00	CM1102 reusable cuff, Child,18-26cm, with connector	Each
	115-020375-00	CM1103 reusable cuff, Adu, 25-35cm, with connector	Each
Intel	115-020376-00	CM1104 reusable cuff, Large Adu,33-47cm,with connector	Each



	connector	
115-020377-00	CM1602 Single-patient cuff liner, Ped, 18-26cm (for CM1102)	Each
115-020378-00	CM1603 Single-patient cuff liner, Adu, 25-35cm (for CM1103)	Each
115-020379-00	CM1604 Single-patient cuff liner, Large Adu, 33-47cm	Each